


CLICKVIEW 7i
**TECHNICAL SPECIFICATIONS –
PACs CONNECTIVITY**


CLICKVIEW 7i

PACS CONNECTIVITY OVERVIEW



CLICKVIEW 7i

PACS ULTRASOUND CONNECTIVITY DETAIL (EXAMPLE) - OPTION 1



ASSUMPTIONS (REQUIRED):


- * Ultrasound Scanner can IMPORT Accession Number or Unique Study ID
- * Ultrasound Scanner can EXPORT Accession Number or Unique Study ID + Ultrasound Data*Link

WORKAROUND:

- * Manual Entry of Accession Number @ Workstation
- * DICOM (See Option #2)

CLICKVIEW 7i

PACS ULTRASOUND CONNECTIVITY DETAIL (EXAMPLE) - OPTION 2



ASSUMPTIONS :

- * CLICKVIEW Server receives Accession Number or USID via DICOM
- * Datalink Transfer of Patient Information and Ultrasound Data

WORKAROUND:

- * Manual Entry of Accession Number @ Workstation